

CELL CULTURE-RELATED PROTOCOLS

C1: Protocol for cell line establishment from embryonic tissues

- 1) Sacrifice mouse using approved procedures
- 2) Immediately swab the animal with 70% alcohol for sterility.
- 3) With sterile forceps and scissors, make a horizontal incision across the abdomen and under the forelegs.
- 4) With sterile scissors and forceps remove the embryos.
- 5) Place embryos in Petri dish containing sterile 1X PBS
- 6) Rinse tissues until the solution is clear.
- 7) Dissect out the appropriate tissues and place them in separate wells containing sterile PBS.
- 8) Under sterile conditions, transfer tissue to separate 25ml sterile flasks containing 5-10ml 0.05% trypsin in HBSS.
- 9) Stir gently in an incubator at 37 degrees for 15 min.
- 10) Using sterile 5ml pipette, pipette up and down the solution to further break-down the remaining tissue.
- 11) Remove the solution containing suspended cells to a sterile 50ml conical plastic centrifuge tube containing 1ml of calf serum per 10ml of suspension to inactivate the trypsin.
- 12) Centrifuge the pooled cells suspension at 1200 rpm for 5 min. and discard the supernatant
- 13) Resuspend the cells pellets in fresh growth media and centrifuge again as in step 1
- 14) Wash the cells several times with growth media until the supernatant is clear.
- 15) Resuspend the final pellets in 5-10ml DMEM containing 10% FCS + antibiotics + high glucose + L-glutamine. Count cells in hemacytometer
- 16) Plate the cells at 10,000,000 cells per 75 cm² flask containing culture media and incubate at 37 degrees in a humidified CO₂ incubator
- 17) Change media after 24 hrs and the every other day until cells reach confluence.

C2: Harvesting Cells with Trypsin/EDTA

- 1) Aspirate media and rinse cells with 5ml sterile 1X PBS.
- 2) Add 4 ml of Tryp/EDTA(0.05% tryp/EDTA) and incubate at 37 degrees for 3-5 min.
- 3) Add equal amount of growth media (RPMI or DMEM) containing 10% FBS-Amp/Strep
- 4) Transfer cells to a sterile tube and centrifuge at 1200 RPM for 5 min. at 4 deg.
- 5) Aspirate media off and add 5-10 ml of growth media (RPMI or DMEM) containing 10%FBS-amp/strep.
- 6) Split appropriately into 75cm² flasks containing 10ml growth media (RPMI or DMEM) containing 10%FBS-amp/strep
- 7) Place cells in the CO₂ incubator.